


Instructions for the use of

idc waste disposal safes

'for contaminated medical materials'



idc Clinical Waste Disposal Safes are in use throughout Australia, New Zealand and the South Pacific Islands.

Thank you . . . for using IDC Disposal Safes!

You have chosen one of the best purpose designed medical waste collection containers available in Australia. You also save storage space because of the 'nesting' design.

How to . . .

Every DISPOSAL SAFE (ill. E) is supplied with one LID (ill. C), one FINGERGUARD (ill. B) one LOCKING RING (ill. D) and one FINAL CLOSURE (ill. A)

Please see reverse for illustrations

Firstly fit the Fingerguard (B) into the Lid (C) opening. Next, screw the Locking Ring (D) onto the thread of the Fingerguard (B). The Fingerguard is now securely 'locked' onto the Lid (C).

The Lid (C) can now be fitted to the Disposal Safe (E - Base Unit). Engage the Lid by pressing down firmly at one point on the Base Unit (E) & then working around the circumference until the Lid is completely secure.

Final Closure (seal when ready for disposal)

The Final Closure (A) 'click fits ' into the Fingerguard (B) apperture. When your Disposal Safe is filled, or ready for disposal simply click the Final Closure into place (Directly into the Lid or Fingerguard, depending). Fitting the Final Closure readies the container for proper disposal according to regional regulations and requirements. Once the Final Closure is sealed into the Fingerguard the Sharps Container cannot be reopened.

Always install and use Fingerguards when the Disposal Safes are used for the disposal of sharps.

idc medical

an Australian Business

Phone: 02 4321 0299 Fax: 02 4321 0288

Catalogue pages and pricelist on: www.idcmedical.com

For your personal protection and the safety of others using the same container . . .

Please remember that:

Polymer products, although puncture resistant, are not entirely puncture proof.

Always follow the guidelines set down by the Infection Control Association. Never fill a container above the indicator fill line. Never force an article into a Disposal Safe. Never, ever deposit syringes with the needle pointing upwards.