

Advantages - Quick Safe

- Secure Fitting Lid:** The Quick Safe base has a patented “tongue garland” closing-rim, & because of this the yellow lid can be fitted easily, quickly & safely. When the lid is fitted onto the base an audible click is heard, signifying a secure watertight closure. Once fitted the lid cannot be removed.
- Two Lid System:** Larger yellow lid fits to base unit. Smaller red lid, which allows daily-use of the sharps container.
- ◆ For daily use the container can be opened & closed as often as required. To temporarily close the Quick Safe, lightly press the front of the day lid into the top opening.
 - ◆ When the container is filled & ready for final closure, firmly press the day lid into the top opening - press down hard around entire perimeter of the day lid until a second ‘click’ is audible. The Quick Safe is now secure & safe for disposal.
- Pull-off Plate:** The pull-off plate (in top opening) enables the secure and problem free screw-off or pull-off of all common cannulas and vacutainer-systems. There is also a star-shaped pull-off device for all common insulin pens.
- Stick-Holder:** Quick Safe containers have formed recess in the bottom of the base. This recess fits the stick-holder bracket (code: QSlock). The Quick Safe can be firmly fixed to any smooth level surface using the Stick holder bracket - just screw the container onto the bracket. In conjunction with this bracket the QS 0.5 container can be fitted onto common syringe –platters/trays.
- Material:** Quick Safe containers are made of environmentally-friendly polypropylene. The red & yellow colours of the containers are free of cadmium.
- Emissions:** Quick Safe containers incorporate no harmful chemicals, such as cadmium or similar. Safe to incinerate with a ‘clean’ burn.



QS 0.5 litre



QS 1.0 litre