

Vario Safe 2.3 Litre Sharps Container

Technical Data Sheet (TDS)

QSvo 2.3

PRODUCT OVERVIEW

The Vario Safe 2.3 Litre Sharps Container is a yellow medical waste sharps container with a variable lid opening designed to accommodate bulky sharps and clinical waste that may not suit smaller apertures. Constructed from polypropylene and supplied as a high-visibility yellow medical waste container, it provides an identifiable, controlled disposal point that supports safe handling of contaminated sharps and clinical waste.



WHY THIS PRODUCT EXISTS

Sharps disposal points can present significant operational and compliance risks, including:

- Needle-stick injuries
- Tampering or unauthorised access to contaminated sharps
- Cross-contamination risks
- Incorrect disposal of sharp or contaminated clinical waste
- Compliance failures where sharps are not controlled at point of use

This container is designed to mitigate these risks by providing a purpose-built sharps disposal point matched to the application, waste type and required capacity.

PRODUCT FUNCTION

The QSvo 2.3 provides controlled containment for sharps and contaminated clinical waste. It supports:

- Variable lid opening for bulky sharps waste
- Configurable deposit aperture by use of an inset
- Water-tight tongue-garland lid-to-base connection
- Clear identification of biohazard and medical waste disposal points
- Compatibility with IDC Medical QS Series sharps waste systems

KEY BENEFITS

- Reduces risk of needle-stick injuries for staff, patients and the public
- Supports safe segregation of contaminated sharps and clinical waste at point of use
- Improves compliance with clinical waste handling procedures
- High-visibility yellow finish increases hazard awareness
- Durable polypropylene construction supports reliable routine use
- Final closure features help prevent access after disposal
- Supports consistent and standardised sharps waste handling procedures

SYSTEM COMPATIBILITY

- Suitable for bulky waste and syringe platter workflows
- Integrates into IDC Medical QS Series sharps disposal infrastructure
- Suitable for clinical waste handling procedures where the nominated capacity and aperture match the waste stream

Vario Safe 2.3 Litre Sharps Container

Technical Data Sheet (TDS)

QSvo 2.3

TECHNICAL SPECIFICATIONS

Specification	Details
Product Type	Seamless Injection Moulded
Product Description	Yellow Medical Waste Vario Safe
Re-order Number	QSvo 2.3
Capacity to Fill-line	2.0 litre
Capacity to Spill-line	2.3 litre
Lid Aperture / Top Opening	95 mm x 65 mm (oval configuration)
Final Closure	Tongue-garland lid-to-base connection
Application	Disposal of bulky sharps waste, e.g. scalpels, used cannulas, disposable syringes and broken glass
Shape	Rectangular / oval profile
Material	Polypropylene
Construction	Seamless injection moulded container with variable lid opening and non-return flap option
Dimensions (Unit)	160 mm (H) x 244 mm (W) x 106 mm (D)
Indentification	Yellow medical waste / sharps container presentation

PACKAGING & LOGISTICS

Specification	Details
Standard Pack Size	30 units per carton
Pack Weight	6 kg
Pack Dimensions	625 mm x 265 mm x 515 mm

Disclaimer:

IDC Medical do not warrant or represent in any way whatsoever that IDC products (Products) meet current standards (including AS23907-2:2019). You expressly acknowledge and agree that you are responsible for ensuring that your use of the Products are compliant for your purposes. You expressly indemnify IDC Medical from and against any claims, actions and all other liability with respect to your use of the Products.

V 202605 - The images and /or information provided in this document are for general guidance only and may be subject to technical changes. We reserve the right to modify product specifications, design, or information where necessary. Images are for illustrative purposes only and may differ from the actual product.